Manual de Usuario / User's Manual / Manual do usuário

DGC-II

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PRECAUTIONS

- The DGC-II should only be used by professionals qualified to use this type of instruments
 - Read this USER'S MANUAL carefully before using the DGC-II for the first time
 - Do not connect any accessory not described in this manual
 - Do not try to open the equipment nor dismantle the probe, as they may suffer irreversible damages and you may receive an electric discharge
 - The detector probe must be treated very gently. Any blow or vibration may cause irreversible damage
 - Before plugging or unplugging the probe from the control unit, make sure that this unit is turned off.
-

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More PRECAUTIONS

- For surgical use, the detector probe must be put into a sterile polyethylene sleeve, with a sterile latex finger cot at its end to avoid direct contact with tissues. In case of possible contamination, it can be cleaned with a cloth dampened in Povidone-Iodine, but not submerged in it. The probe can also be sterilized with ethylene-oxide.
 - Do not attempt to use the DGC-II during the battery loading cycle, as it will be blocked.
 - Do not use the DGC-II for continuous periods of over 60 minutes.
 - Do not wet the equipment nor handle it with wet hands.
 - Handling of radioactive material must be done following the regulations of the Nuclear Regulatory Authority (or the corresponding authority in each country)
- Do not use the DGC-II for any purpose other than that for which it was designed.
-

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Introduction

DGC-II is a system specifically designed to identify high captation *hot spots* of radioactive material during surgery.

DGC-II is composed of a detector probe, a visualization and control unit and a cable that interconnects both.

Your detector probe has enough space resolution and count effectiveness needed to locate small *hot spots* such as Sentinel Nodes or other small *hot spots* sought in radio-guided surgical procedures.

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The control and data visualization unit receives signals from the detector probe, analyzes them and shows, in digital form, the count rate by means of a numerical display and analogically via an arrangement of LEDs and an audible signal, whose beep frequency is proportional to the count rate. This unit works with an internal, rechargeable battery and should therefore not be plugged into the electric network to operate. To recharge the batteries, connect the control unit to an external power source (not provided by the manufacturer), that can supply 18 VAC @ 1 A. The operation of the instrument is blocked while it is plugged to the external source.

DGC-II can detect gamma radiation from different isotopes and is pre-set by manufacturer to detect Tc-99m, as it is the most widely used for this application.

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 - Accessories not supplied by the manufacturer:
 - External source for supply of 18 VAC @ 1 A.
 - Carrier bag.
-

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Description of the system

Detector probe and interface cable

The detector probe is contained in a stainless steel housing; the detector head is on one end and the interface cable connector on the other.

The head contains a CdZnTe crystal of 5 x 5 x 3 mm placed inside a shielding- collimator tungsten assembly. The central shaft of the probe contains a low-noise pre-amplifier that arranges the signals that are produced in the crystal when the gamma rays from the radioactive material act on it.

The interface cable has two connectors. The metal one is connected to the detector probe and the plastic one (gray) is connected to the control and visualization unit.

To connect the interface cable to the detector probe, both red marks must meet (the one on the metal connector of the cable with the corresponding one of the probe) and the connector careful inserted. The connectors allow only one position and their connection must never be forced. To prolong the life of the cable and the connectors, they should be carefully rolled for storage, avoiding bends and use of force.

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Control and visualization unit

The control and data visualization unit receives signals from the detector probe, analyzes them and shows, in digital form, the count rate by means of a numerical display and analogically via an arrangement of LEDs and an audible signal, whose beep frequency or tone is proportional to the count rate. This unit works with a rechargeable internal battery and should therefore not be plugged into the electric network to operate. To recharge the batteries, connect the control unit to an external power source (not provided by the manufacturer), that can supply 18 VAC @ 1 A. The operation of the instrument is blocked while it is plugged to the external source.

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The following controls and/or indicators are on the front panel:

Power On key

State of battery indicator

Probe connector
Scale range Selector
Numeric display
Analog display (Luminous bars)
Integration time selector
Window selector
Audio mode selection key

The following are on the rear panel:

Connector for battery charger
Volume adjustment knob

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Power on key

It has two positions “0” and “1”

In position “0” the equipment is off. To recharge the batteries, the power on key must be in this position.

In position “1” the equipment is on. When the key is turned to this position, the **numeric display** will light up.

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State of battery indicator

This luminous indicator has different states:

With the key in position **OFF** and the external power source connected:

- ◇ **GREEN - STILL:** Charging battery
- ◇ **GREEN - BLINKING:** Charge completed
- ◇ **OFF:** External power source faulty or disconnected

With the key in position **ON** and the external power source disconnected:

- ◇ **RED:** Low battery. Indicates batteries need to be recharged. The equipment has autonomy for half an hour after this indicator lights up.
 - ◇ **OFF:** Normal operation
-

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The battery of the DGC-II has a life span of 2 years. After that time, please contact the factory for its replacement. Due to its nature, the battery of the DGC-II does not accumulate memory, so the user can recharge it as often as needed. Although the DGC-II has an internal battery charge control, the system should not be charged for periods of more than 16 hours. If after that time the full charge has not been achieved, have the manufacturer check the equipment.

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Connector for the probe

The plastic connector of the detector probe cable is inserted in this connector.

For this connection, insert the plastic connector of the cable with the arrow pointing upwards on the control and visualization unit. A CLICK must be heard, to make sure that the connector is properly in place. The connectors allow only one position and must never be forced. To disconnect the probe, slide the broad part of the connector backwards. This will unlock it and make it easy to remove.

BEFORE CONNECTING OR DISCONNECTING THE PROBE, MAKE SURE THAT THE INSTRUMENT IS TURNED OFF.

DO NOT ATTEMPT TO REMOVE THE CONNECTOR BY PULLING ON THE CABLE OR UNSCREWING, AS THIS WILL DAMAGE IT IRREVERSIBLY.

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Scale range Selector

This knob allows selection of the different scale ranges for the ANALOG DISPLAY (luminous bars) and the audible signal (series of beeps). The available scale ranges are 200; 500; 2000; 5000 and 20000 cps.

The position of this control does not affect the data shown in the NUMERIC DISPLAY. It only allows adjustment of the luminous bars and the beeps to the range of the job. Always begin with the lowest range (200) and change to the immediately higher one should there be saturation in that range (all the luminous bars remain lit).

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NUMERIC DISPLAY

In this display the control and visualization unit shows the count rate in digital form, expressed in **cps** (counts per second).

The maximum reading rate is 20000 cps. Beyond this, the system will be saturated and will indicate this by showing all the digits blinking.

Should it be necessary to work with high count rates, and it were noticed that the instrument is saturated, set the WINDOW key in NARROW position and try again.

Analog display (Luminous bars)

This analog display allows a rapid view of the variations of the count rate. Its sensitivity is adjusted with the SCALE RANGE selector (see SCALE RANGE SELECTOR).

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Integration time selector

The integration time is the period during which the control unit accumulates. Taking into account that the radioactive disintegration phenomenon is a random one, short integration periods will allow fast sweepings with the probe, but low count rates will result in a fluctuating reading.

Three different integration times can be selected:

“1”: Short (suitable for fast sweeps)

“2”: Intermediate (recommended)

“3”: Slow (assures greater reading stability; use only to measure fixed points)

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Window selector

The measuring window represents the part of the gamma spectrum that will be taken into account during the analysis of the pulses detected by the probe. For this purpose, two pulse discriminators are set (power thresholds): a low and a high one. Only those pulses that pass through the preset window will be taken into account in the analysis.

DGC-II is delivered by manufacturer adjusted for Tc-99m measuring.

The **wide** position (key down) captures pulses between 100 and 180 keV and is the one **recommended for use with Tc-99m**.

The **narrow** position (key up) captures pulses between 130 and 150 keV and its use is only recommended for very high count rates to avoid saturation of the DGC-II.

The **open** position sets the lower discriminator at 20keV and ignores the higher one. Bear in mind that in this position, spurious pulses can affect the spatial resolution.

If you want to **use I-125**, you **must** select the **open** position.

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Audio mode selector key

Pressing this key down, the audio will be silenced.

Volume adjustment knob

Audio volume can be adjusted with this knob.

Connector for battery charger

This is the power source entrance for the internal charger of the rechargeable batteries.

To recharge the batteries, connect the control unit to an external charger (not provided by the manufacturer) that can supply 18 VAC @ 1 A. While the instrument is connected to the external charger, its operation is blocked.

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OPERATION

General Considerations

- DGC-II is delivered calibrated for **Tc-99m** measuring.
 - The system may only be operated by qualified personnel for this type of work.
 - Handling of radioactive material must be done following the regulations of the Nuclear Regulatory Authority (or the corresponding authority in each country)
 - To be used, the detector probe must be put into a sterile polyethylene sleeve, with a sterile latex finger cot at its end to avoid direct contact with tissues potentially contaminated with radioactive matter. Should there be any contamination, the probe can be cleaned with a cloth dampened in Povidone-Iodine, but not submerged in it.
-

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Connecting the probe.

Make sure that the equipment is turned off.

Connect the detector probe to the control and visualization unit, introducing the plastic connector of the probe to the one labeled PROBE

A click will indicate that the probe is properly connected.

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Adjusting the SCALE RANGE

Put the SCALE RANGE SELECTOR in position 200 when starting a sweep. If you notice that the equipment becomes saturated in that position, move the SCALE RANGE selector to the next immediate level until reaching the required sensitivity.

SELECTING THE WINDOWS

Put the WINDOW SELECTOR in position WIDE, pressing down the key. Only in case you notice that the equipment is saturated due to a high count rate, change to the NARROW position.

The OPEN WINDOW position should only be used in case of very low count rates, or when you're using I-125 as tracer, as in this position spurious pulses can affect the measurement.

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SELECTING INTEGRATION TIME

Set the INTEGRATION TIME SELECTOR in position “2” (INTERMEDIATE). Only in case of needing to measure a fixed point with greater reading stability should the selector be put in position “3”. Position “1” is reserved for very fast sweeps.

Adjusting the audio

Use the AUDIO MODE selector key and the volume control on the rear panel to adjust audio and volume.

The DGC-II is ready to measure. Should you observe that the equipment is saturated, pass the SCALE RANGE control to the immediately higher position.

Important Note: To be used, the detector probe must be put into a sterile polyethylene sleeve, with a similarly sterile latex finger cot at its end to avoid direct contact with the tissues. Should there be any contamination, the probe can be cleaned with a cloth dampened in Povidone-Iodine, but not submerged in it.

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After using the DGC-II

Proceed as follows:

- Turn the equipment off, putting the power on key in position "0" .
- Disconnect the probe from the control and visualization unit. To disconnect the probe, slide the broad part of the connector backwards. This will unlock it and make it easy to remove.

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Wipe the probe and its cable and the cabinet of the visualization and control unit with a cloth soaked in alcohol 70°.

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Recharging the battery

This unit works with a rechargeable internal battery and should therefore not be plugged into the electric network to operate.

To recharge the battery, connect the control unit to an external power source (not provided by the manufacturer), that can supply 18 VAC @ 1 A. The operation of the instrument is blocked while it is plugged to the external source.

The battery of the DGC-II has a life span of 2 years. After that time, please contact the factory for its replacement. Due to its nature, the battery of the DGC-II does not accumulate memory, so the user can recharge it as often as needed. Although the DGC-II has an internal battery charge control, the system should not be left charging for over 16 hours. If after that time the full charge has not been achieved, have the manufacturer check the equipment.

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CHECKING THE EFFICIENCY OF THE SYSTEM

As in the applications of the DGC-II it is not necessary to know the absolute activity of radioactivity in the points of interest, but only to detect high pick-up areas, THE DGC-II SYSTEM DOES NOT REQUIRE CALIBRATION for its use for the purpose for which it was designed.

However, we recommend making periodic CHECKS OF THE SYSTEM'S EFFICIENCY at least every six months.

Required Materials

Sealed Co-57 Source with an activity between 50 and 100 μ Curies
Scale in cm.

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Procedure

Adjust the parameters of the DGC-II in the following position:

WINDOW: WIDE

INTEGRATION TIME: "5"

Position the active face of the probe towards the Co-57 source, at a distance of $130 \text{ mm} \pm 5 \text{ mm}$ until you obtain the maximum reading.

Under these conditions, there should be a reading between 2 and 3 cps/ μ Curie

Should the reading not be within that range, contact the manufacturer.

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TECHNICAL SPECIFICATIONS

Control and visualization Unit

- Power: Internal rechargeable batteries
- Battery charger powered by an external source (not provided by the manufacturer) that can supply 18 VAC @ 1 A
- Maximum count rate: 20000 cps (Counts per second)
- Temperature range during operation: + 15°C a + 40°C
- Storage temperature range: + 5°C to + 40°C
- Relative humidity: less than 80% (non-condensing)
- Preset isotope: Tc-99m

Nominal autonomy of internal battery: > 7 hours with maximum charge.

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- Height: 12,5 cm
- Width: 28,5 cm
- Depth: 30 cm
- Weight: 4 kg

Detector Probe

- Type of Detector: Semiconductor of CdTe (Zn)
 - Dimensions of the detector: 5 x 5 x 3 mm
 - Low noise pre-amplifier built into central shaft
 - Energy range: 20 to 200 keV
 - Type of Collimator: Internal, made of tungsten
-

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- Storage temperature range: + 5°C to + 40°C
- Housing material: Stainless steel
- Dimensions (without cable):
 - Total length: 160 mm
 - Diameter of the head: 13 mm
 - Central shaft diameter: 10 mm
- Angle of the head, in relation to central shaft: 35° (simplifies access in small incisions)
- Weight: 70 g
- Relative humidity < 80% (not condensing)

All dimensions and parameters stated in this USER'S MANUAL are nominal and may be changed by the manufacturer without prior notice.

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